3130170

CE 0197

CARE AND USE OF YOUR PEAK FLOW METER

It is important to keep your Peak Flow Meter clean when in use. This should be cleaned at least once every six months and thoroughly if it has not been used for some time, fungal growths may otherwise occur. The Mouthpiece is detachable and should be disinfected before each use using a proprietary disinfectant solution such as Dettol. The OMRON PFM20 Peak Flow Meter is designed for single patient use.

Technical Information The OMRON PFM20 Peak Flow Meter conforms to ISO 23747 and the NAEPP Guidelines. Accuracy: Better than 10% or 10l/min. Readings expressed as BTPS values. Maximum resistance to flow 0.35 kPa/l/s. Store in a dry place.

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YOUR PEAK FLOW METER Your OMRON PFM20 Peak Flow Meter will help you to monitor your asthma by measuring your peak flow. This leaflet explains how to use your Peak Flow Meter, and how to record the readings. With careful use the OMRON PFM20 would be expected to last for two years.

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



1 If required insert a mouthpiece into the meter. Ensure the pointer is set at Zero.


2 Hold the Peak Flow Meter so that your fingers are clear of the scale and slot. Do not obstruct the holes at the end of the Peak Flow Meter, or the vents at the side of the mouthpiece.

3 Stand up if possible. Take a deep breath, place the Peak Flow Meter in the mouth and hold horizontally, closing the lips around the mouthpiece, then blow as hard and as fast as you can. Coughing and spitting into your meter will adversely affect your reading and should be avoided.

4 Note the number on the scale indicated by the pointer. If unusual readings are obtained and the meter does not appear broken - seek advice from your healthcare provider.

5 Return the pointer to zero and repeat the procedure twice more to obtain three readings. Mark the highest of the three readings on the chart as shown overleaf.

CLEANING YOUR PEAK FLOW METER

Immerse Peak Flow Meter in warm (but not hot) mild detergent solution for 2-3 minutes (maximum 5 minutes). Agitate the meter to ensure thorough cleaning.

Rinse in clean warm water and shake gently to remove any excess water. Allow to dry thoroughly before using again.

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
**16 HOW TO RECORD YOUR PEAK FLOW**

It is important to keep a record of your peak flow. Readings should be taken twice a day, first thing in the morning, and at bed time, or as directed by your physician. Use the chart provided to record the highest reading obtained each time. Mark your reading on the chart with a dot corresponding to the highest reading obtained for each of the morning, and bed time readings every day. By joining the dots (see diagram) you can see how your peak flow varies from day to day. You should also measure and note your peak flow if you suddenly become wheezy or breathless at any time of the day or night.

The use of a 3 zone monitoring system set up like traffic lights uses coloured plastic clips (green, yellow, red) and has proven useful to patients watching their asthma with peak flow readings. Once your green, yellow and red zones have been determined, record their values in the space above. Then set the zone markers on your meter to correspond with those values. Space has also been provided above to record the actions to be taken, from your treatment plan, based on your peak flow readings.

The general principles of the zone monitoring system are as follows, and would normally be set by your Doctor or Nurse:

Zone	Condition	Reading (% of best)
Green Zone	Go ahead with normal activities	80-100% of best reading
Yellow Zone	Caution	60-80% of best reading
Red Zone	STOP	Less than 50% of best reading

*Example: If you and your physician or other licensed healthcare professional determine your individual best reading is 500 LPM, then your Green Zone is 400-500 LPM (80-100%). Set the green marker at 500. Your Yellow Zone is 250-400 LPM (50-80%). Set the yellow marker at 400. Your Red Zone is 0-250 (less than 50%). Set the red marker at 250.

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