MyColor

DIGITAL COMPACT THERMOMETER

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We congratulate you on your purchase of this thermometer. To correctly utilise all the advantages of this product, you should familiarise yourself with these operating instructions before using the thermometer the first time. Then keep these instructions in a safe place for future reference.

SPECIAL FEATURES

- **Maximum thermometer:** in each case the highest measured temperature is displayed until the appliance switches off.
- **Automatic self-test:** guarantees a specified measuring accuracy for each temperature reading! «ERR» (error) is immediately indicated if any deviation is registered.
- **Temperature measurement** within 60-90 seconds depending on measuring method used.
- **Digital display** for temperature changes as minor as 0,1°C
- **Memory:** automatic storage of last measured temperature.
- **Automatic switchoff:** after approx. 10 minutes – for longer battery life.
- **100 % waterproof:** can be completely immersed for cleaning.
- **Battery:** Silver oxide high drain. (Non disposable)
- **Measurement in:** Celsius.

IMPORTANT SAFETY INSTRUCTIONS

- The thermometer is suitable exclusively for measurement of the body temperature!
- The battery cover is sonic-welded with the main housing of the thermometer; the battery is well sealed inside the thermometer to reduce the risk to be swallowed by children.
- Protect the thermometer from hitting and striking!
- Avoid ambient temperatures above 60°C. NEVER boil the thermometer!
- Use only the commercial disinfectants listed in the section «cleaning and disinfecting» to clean the thermometer. The thermometer must be intact when immersed in liquid disinfectant.
- **Warning:** The thermometer contains a long-life battery. For disposal, please comply with local regulations for special waste.
- The manufacturer recommends to verify the accuracy by an authorised laboratory every 2 years.

TURNING ON THE THERMOMETER

To turn on the thermometer, push the button on top; the display test is performed. Then at an ambient temperature of less than 32°C, an «L» and a flashing «°C» appear at the upper right of the display. The thermometer is now ready for use.

FUNCTION TEST

Proper functioning of the thermometer is tested automatically each time it is turned on. If a malfunction is detected (measurement inaccuracy), this is indicated by «ERR» on the display, and a measurement becomes impossible. In this case, the thermometer must be replaced.

USING THE THERMOMETER

During a measurement, the current temperature is displayed continuously. The «°C» symbol flashes during this period of time. If the «°C» is no longer flashing, this means that the measured increase in temperature is now less than 0,1 °C in 16 seconds and that the thermometer is ready for reading. To prolong the life of the battery, turn it off after use by pressing the operating button briefly. Otherwise the thermometer will automatically turn off after about 10 minutes.

STORAGE OF MEASURED VALUES

If you hold the button down for more than 3 seconds when turning the thermometer on, the maximum temperature stored automatically at the last measurement will be displayed. At the same time, a small «M» for memory will appear at the right of the display. After 3 more seconds, this value will disappear and the thermometer will switch to the normal measurement mode. Then the stored value is erased.

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MEASURING METHODS

In the mouth (oral)

Position the thermometer in one of the two pockets under the tongue, to the left or right of the root of the tongue. The measuring sensor must be in good contact with the tissue. Close your mouth and breath evenly through the nose to prevent the measurement from being influenced by inhaled/exhaled air.
Approx. measuring time: 1 minute.

In the anus (rectal)

This is the most reliable measuring method, and is especially appropriate for infants and small children. Carefully insert the tip of the thermometer 2 to 3 cm into the anal aperture.
Approx. measuring time: 1 minute!

In the armpit (axillary)

From a medical viewpoint, this method will always provide inaccurate readings, and should not be used if precise measurements are required.
Approx. measuring time: 5 to 10 minutes!

CLEANING AND DISINFECTING

The thermometer is 100% waterproof and may be immersed in water or disinfectant solution for cleaning with no problem. Then wipe it off with a moist towel. Please use one of the following disinfectants for disinfecting, following the manufacturer's instructions in each case.

Name of disinfectant	Manufacturer	Dilution with water	Max. duration Method
Isopropyl alcohol 70%		95,0%	Immerse for 24 hours
Gigasept FF	Schilke & Mayr GmbH	5,0%	Immerse for 24 hours
Lysoformin 3000	Rosemann GmbH	6,0%	Immerse for 24 hours
Aseptisol	Bode Chemie	4,0%	Immerse for 24 hours

Note: Any other cleaning solution, dilution or cleaning method can cause a malfunction of the thermometer!
The thermometer must never be boiled!

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WARRANTY

We grant you a lifetime warranty after the date of purchase. Any damage caused by improper handling shall not be covered by the warranty. Batteries and packaging are also excluded from the warranty. All other damage claims excluded. A warranty claim must be submitted within the warranty period. Please package your defective thermometer well and send with sufficient postage to:

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